

DIRECTION FOR USE: To clean, rinse, disinfect and remove protein daily, these “rub and rinse” steps must be completed:

Preparations:

Do not use if sealed dropper tip is pierced or damaged at time of purchase.
To open bottle, remove cap and protective ring. Replace cap on bottle. Turn cap (containing sharp pin) clockwise to pierce bottle’s sealed dropper tip.
Thoroughly wash and clean hands before handling contact lenses

INSTRUCTIONS FOR USE:

(Rub and Rinse Directions are for all contact lens replacement schedules)

Step 1: Cleaning and Protein Removal. Apply 2 to 3 drops of PrimaLens® Multi-Purpose Solution on each side of your contact lens. Gently rub lens from the center to the edges in the palm of your hand using the forefinger of the other hand for approximately 20 seconds. Be sure to clean both sides. No separate daily cleaner necessary.

Step 2: Rinsing. Thoroughly rinse each side of each lens for 10 seconds with PrimaLens® Multi-Purpose Solution to remove all debris.

Step 3: Disinfection and Storage. Fill your contact lens case with PrimaLens® Multi-Purpose Solution and place lenses in the correct side. Close case and tighten. Allow lenses to soak at least 4 hours.

Your lenses are now ready to wear. No saline rinse is necessary. If any debris remains on contact lenses, rinse with PrimaLens® Multi-Purpose Contact Lens Solution. Always follow your Eye Care professional’s instructions. Based on your individual tear chemistry and lens-wearing schedule, your Eye Care professional may recommend the use of a separate enzymatic cleaner before step 1.

Storage: Using a clean case, completely submerge your lenses in PrimaLens® Multi-Purpose Contact Lens Solution to store them. Lenses may be stored in the unopened case until ready to wear, up to a maximum of 30 days. If you store your lenses for longer periods of time, they must be cleaned and disinfected with fresh PrimaLens® Multi-Purpose Contact Lens Solution every 30 days prior to wearing.

Protein Removal: PrimaLens® Multi-Purpose Contact Lens Solution removes protein from your contact lenses while you clean and disinfect your lenses daily.

Shelf life: Please discard bottle 60 days after opening.

Can be used for all Soft Contact Lenses to:

Clean • Rinse • Store • Disinfect • Remove Protein Deposits

No saline solution is required with the PrimaLens® Multi-Purpose Contact Lens Solution.

NOTE: Always use fresh PrimaLens® Multi-Purpose Contact Lens Solution & never reuse the solution.

WARNINGS: To avoid contamination, nothing should ever touch the dropper tip. Keep container tightly closed when not in use.

CONTAINS: PrimaLens® Multi-Purpose Contact Lens Solution is a sterile, isotonic, buffered, solution containing Boric Acid, Sodium Chloride, Hydroxypropyl Methylcellulose (HPMC) as a lubricant, Poloxamer 407 as a surfactant, Disodium Edetate as chelating agent, Purified Water and preserved with Polyhexanide.

If you are allergic to any of the above ingredients do not use this product

TAMPER EVIDENT: For your protection, this bottle has a sealed dropper tip. Do not use if sealed dropper tip is pierced or damaged at time of purchase.

Keep Out of Reach of Children

Distributed by:

Molecular Health Commerce, Inc.

San Diego, CA 92122

For further information contact

1-866-766-MHCI or go to PrimaLensMPS.com

Made in Germany